

Collagen I Mouse Monoclonal Antibody(6E11)

Catalog No: #44271

Package Size: #44271-1 50ul #44271-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Collagen I Mouse Monoclonal Antibody(6E11)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	B'B 6E11
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Collagen I protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of Collagen I
Target Name	Collagen I
Other Names	Collagen alpha-1(I) chain (Alpha-1 type I collagen)
Accession No.	Swiss-Prot:P02452Gene ID:1277
Uniprot	P02452
GeneID	1277;
Calculated MW	139kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Collagen I Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Skin Carcinoma Tissue using Collagen I Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only