

Cytochrome P450 2C19

Catalog No: #AP82092



Package Size: #AP82092-1 50ug #AP82092-2 100ug #AP82092-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 2C19
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	(R)-limonene 6-monooxygenase,(S)-limonene 6-monooxygenase,(S)-limonene 7-monooxygenase,CYPIIC17,CYPIIC19,Cytochrome P450-11A,Cytochrome P450-254C,Mephenytoin 4-hydroxylase
Accession No.	Uniprot ID:P33261
Uniprot	P33261
GeneID	1557;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CYP2C19

Note: This product is for in vitro research use only