

Cysteine-rich secretory protein LCCL domain-containing 2



Catalog No: #AP82093

Orders: order@signalwayantibody.com

Package Size: #AP82093-1 50ug #AP82093-2 100ug #AP82093-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Cysteine-rich secretory protein LCCL domain-containing 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:A6QLZ7
Uniprot	A6QLZ7
GeneID	505329;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CRISPLD2

Note: This product is for in vitro research use only