

Casein kinase II subunit beta

Catalog No: #AP82096



Package Size: #AP82096-1 50ug #AP82096-2 100ug #AP82096-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Casein kinase II subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CK II beta, Phosvitin, Protein G5a, CK2N, G5A
Accession No.	Uniprot ID:P67870
Uniprot	P67870
GeneID	1460;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CSNK2B

Note: This product is for in vitro research use only