

Collagen alpha-1(XVII) chain

Catalog No: #AP82098



Package Size: #AP82098-1 50ug #AP82098-2 100ug #AP82098-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(XVII) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	180 kDa bullous pemphigoid antigen 2,Bullous pemphigoid antigen 2,BP180,BPAG2
Accession No.	Uniprot ID:A6QPB3
Uniprot	A6QPB3
GeneID	513804;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COL17A1

Note: This product is for in vitro research use only