

Cytochrome c oxidase subunit 3

Catalog No: #AP82099



Package Size: #AP82099-1 50ug #AP82099-2 100ug #AP82099-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome c oxidase subunit 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Cytochrome c oxidase polypeptide III, COIII, COXIII, MTCO3
Accession No.	Uniprot ID:Q9ZZ61
Uniprot	Q9ZZ61
GeneID	804480;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MT-CO3

Note: This product is for in vitro research use only