

Beta-Ala-His dipeptidase

Catalog No: #AP82100



Package Size: #AP82100-1 50ug #AP82100-2 100ug #AP82100-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-Ala-His dipeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CNDP dipeptidase 1,Carnosine dipeptidase 1,Cn1
Accession No.	Uniprot ID:Q8BUG2
Uniprot	Q8BUG2
GeneID	338403;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cndp1

Note: This product is for in vitro research use only