

Complement component C8 beta chain

Catalog No: #AP82101



Package Size: #AP82101-1 50ug #AP82101-2 100ug #AP82101-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement component C8 beta chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Complement component 8 subunit beta
Accession No.	Uniprot ID:Q8BH35
Uniprot	Q8BH35
GeneID	110382;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: C8b

Note: This product is for in vitro research use only