

Macrophage colony-stimulating factor 1

Catalog No: #AP82103



Package Size: #AP82103-1 50ug #AP82103-2 100ug #AP82103-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Macrophage colony-stimulating factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CSF-1,Lanimostim
Accession No.	Uniprot ID:P09603
Uniprot	P09603
GeneID	1435;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CSF1

Note: This product is for in vitro research use only