

Neutrophil defensin 3

Catalog No: #AP82106



Package Size: #AP82106-1 50ug #AP82106-2 100ug #AP82106-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neutrophil defensin 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Defensin, alpha 3,HNP-3,HP-3,DEF3
Accession No.	Uniprot ID:P59666
Uniprot	P59666
GeneID	1668;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DEFA3

Note: This product is for in vitro research use only