

## C-X-C chemokine receptor type 7

Catalog No: #AP82107



Package Size: #AP82107-1 50ug #AP82107-2 100ug #AP82107-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-X-C chemokine receptor type 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	CXC-R7,CXCR-7,Chemokine orphan receptor 1,G-protein coupled receptor RDC1,RDC-1,CMKOR1,RDC1
Accession No.	Uniprot ID:P11613
Uniprot	P11613
GeneID	403964;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CXCR7

Note: This product is for in vitro research use only