

Cytoplasmic FMR1-interacting protein 2

Catalog No: #AP82109



Package Size: #AP82109-1 50ug #AP82109-2 100ug #AP82109-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytoplasmic FMR1-interacting protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	p53-inducible protein 121,Kiaa1168,Pir121
Accession No.	Uniprot ID:Q5SQX6
Uniprot	Q5SQX6
GeneID	76884;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cyfp2

Note: This product is for in vitro research use only