

Beta-defensin 114

Catalog No: #AP82113



Package Size: #AP82113-1 50ug #AP82113-2 100ug #AP82113-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-defensin 114
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Beta-defensin 14,DEFB-14,Defensin, beta 114,DEFB14
Accession No.	Uniprot ID:Q30KQ6
Uniprot	Q30KQ6
GeneID	245928;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DEFB114

Note: This product is for in vitro research use only