

C-reactive protein

Catalog No: #AP82114



Package Size: #AP82114-1 50ug #AP82114-2 100ug #AP82114-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-reactive protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ptx1
Accession No.	Uniprot ID:P14847
Uniprot	P14847
GeneID	12944;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Crp

Note: This product is for in vitro research use only