

Atrial natriuretic peptide-converting enzyme

Catalog No: #AP82116



Package Size: #AP82116-1 50ug #AP82116-2 100ug #AP82116-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Atrial natriuretic peptide-converting enzyme
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Corin,Heart-specific serine proteinase ATC2,Pro-ANP-converting enzyme,Transmembrane protease serine 10,CRN,TMPRSS10
Accession No.	Uniprot ID:Q9Y5Q5
Uniprot	Q9Y5Q5
GeneID	10699;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CORIN

Note: This product is for in vitro research use only