

Cytochrome c oxidase subunit 3

Catalog No: #AP82117



Package Size: #AP82117-1 50ug #AP82117-2 100ug #AP82117-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome c oxidase subunit 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytochrome c oxidase polypeptide III, Cooiii, mt-Co3
Accession No.	Uniprot ID:P05505
Uniprot	P05505
GeneID	26204;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mtco3

Note: This product is for in vitro research use only