

CUB and zona pellucida-like domain-containing protein 1

Catalog No: #AP82118

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP82118-1 50ug #AP82118-2 100ug #AP82118-3 1mg

Description

Product Name	CUB and zona pellucida-like domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CUB and ZP domain-containing protein 1, Transmembrane protein UO-44, UNQ224/PRO257
Accession No.	Uniprot ID:Q86UP6
Uniprot	Q86UP6
GeneID	50624;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CUZD1

Note: This product is for in vitro research use only