

C-X-C chemokine receptor type 4

Catalog No: #AP82119



Package Size: #AP82119-1 50ug #AP82119-2 100ug #AP82119-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C chemokine receptor type 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CXC-R4,Fusin,LCR1,Leukocyte-derived seven transmembrane domain receptor,LESTR,Stromal cell-derived factor 1 receptor,SDF-1 receptor,CD184
Accession No.	Uniprot ID:P25930
Uniprot	P25930
GeneID	281736;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CXCR4

Note: This product is for in vitro research use only