

Sterol 26-hydroxylase, mitochondrial

Catalog No: #AP82122



Package Size: #AP82122-1 50ug #AP82122-2 100ug #AP82122-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sterol 26-hydroxylase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	5-beta-cholestane-3-alpha,7-alpha,12-alpha-triol 27-hydroxylase,Cytochrome P-450C27/25,Cytochrome P450 27,Sterol 27-hydroxylase,Vitamin D(3) 25-hydroxylase,Cyp27
Accession No.	Uniprot ID:P17178
Uniprot	P17178
GeneID	301517;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cyp27a1

Note: This product is for in vitro research use only