

Carnitine O-palmitoyltransferase 1, liver isoform

Catalog No: #AP82123



Package Size: #AP82123-1 50ug #AP82123-2 100ug #AP82123-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carnitine O-palmitoyltransferase 1, liver isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CPT1-L, Carnitine O-palmitoyltransferase I, liver isoform, CPT I, Carnitine palmitoyltransferase 1A, CPT1
Accession No.	Uniprot ID:P50416
Uniprot	P50416
GeneID	1374;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CPT1A

Note: This product is for in vitro research use only