

Corticoliberin

Catalog No: #AP82124



Package Size: #AP82124-1 50ug #AP82124-2 100ug #AP82124-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Corticoliberin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Corticotropin-releasing factor,CRF,Corticotropin-releasing hormone
Accession No.	Uniprot ID:P01143
Uniprot	P01143
GeneID	81648;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Crh

Note: This product is for in vitro research use only