

Chromogranin-A

Catalog No: #AP82126



Package Size: #AP82126-1 50ug #AP82126-2 100ug #AP82126-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chromogranin-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CgA
Accession No.	Uniprot ID:P04404
Uniprot	P04404
GeneID	397540;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CHGA

Note: This product is for in vitro research use only