

## Claudin-3

Catalog No: #AP82134



Package Size: #AP82134-1 50ug #AP82134-2 100ug #AP82134-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Claudin-3  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Hu   |
| Other Names        | Clostridium perfringens enterotoxin receptor 2,CPE-R 2,Rat ventral prostate.1 protein<br>homolog,hRVP1,C7orf1,CPETR2 |
| Accession No.      | Uniprot ID:O15551  |
| Uniprot            | O15551   |
| GeneID             | 1365;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)   |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CLDN3

Note: This product is for in vitro research use only