

Cytochrome c oxidase subunit 5A, mitochondrial

Catalog No: #AP82135



Package Size: #AP82135-1 50ug #AP82135-2 100ug #AP82135-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome c oxidase subunit 5A, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cytochrome c oxidase polypeptide Va
Accession No.	Uniprot ID:P12787
Uniprot	P12787
GeneID	12858;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cox5a

Note: This product is for in vitro research use only