

Coenzyme Q-binding protein COQ10 homolog B, mitochondrial



Catalog No: #AP82137

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AP82137-1 50ug #AP82137-2 100ug #AP82137-3 1mg

Description

Product Name	Coenzyme Q-binding protein COQ10 homolog B, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q3THF9
Uniprot	Q3THF9
GeneID	67876;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Coq10b

Note: This product is for in vitro research use only