

Alpha-defensin 1

Catalog No: #AP82138



Package Size: #AP82138-1 50ug #AP82138-2 100ug #AP82138-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-defensin 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cryptdin-1,Defensin-related cryptdin peptide,Defcr,Defcr1
Accession No.	Uniprot ID:P11477
Uniprot	P11477
GeneID	13216;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Defa1

Note: This product is for in vitro research use only