

Testosterone 17-beta-dehydrogenase 3

Catalog No: #AP82141



Package Size: #AP82141-1 50ug #AP82141-2 100ug #AP82141-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Testosterone 17-beta-dehydrogenase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	17-beta-hydroxysteroid dehydrogenase type 3, 17-beta-HSD 3, Testicular 17-beta-hydroxysteroid dehydrogenase, Edh17b3
Accession No.	Uniprot ID:O54939
Uniprot	O54939
GeneID	117182;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hsd17b3

Note: This product is for in vitro research use only