

Cytochrome P450 1A2

Catalog No: #AP82142



Package Size: #AP82142-1 50ug #AP82142-2 100ug #AP82142-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 1A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	CYP1A2,Cytochrome P450-D2,DAH2
Accession No.	Uniprot ID:P56592
Uniprot	P56592
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CYP1A2

Note: This product is for in vitro research use only