

C-X-C chemokine receptor type 3

Catalog No: #AP82143



Package Size: #AP82143-1 50ug #AP82143-2 100ug #AP82143-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C chemokine receptor type 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CXC-R3,CKR-L2,G protein-coupled receptor 9,Interferon-inducible protein 10 receptor,IP-10 receptor,CD183,GPR9
Accession No.	Uniprot ID:P49682
Uniprot	P49682
GeneID	2833;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CXCR3

Note: This product is for in vitro research use only