

## Carnitine O-palmitoyltransferase 1, liver isoform

Catalog No: #AP82146



Package Size: #AP82146-1 50ug #AP82146-2 100ug #AP82146-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Carnitine O-palmitoyltransferase 1, liver isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CPT1-L, Carnitine O-palmitoyltransferase I, liver isoform, CPT I, Carnitine palmitoyltransferase 1A, Cpt-1, Cpt1
Accession No.	Uniprot ID:P32198
Uniprot	P32198
GeneID	25757;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Cpt1a

Note: This product is for in vitro research use only