

Cerebral dopamine neurotrophic factor

Catalog No: #AP82148



Package Size: #AP82148-1 50ug #AP82148-2 100ug #AP82148-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cerebral dopamine neurotrophic factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ARMET-like protein 1, Conserved dopamine neurotrophic factor, Armet1
Accession No.	Uniprot ID:Q8CC36
Uniprot	Q8CC36
GeneID	227526;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cdnf

Note: This product is for in vitro research use only