

C-X-C chemokine receptor type 4

Catalog No: #AP82156

Package Size: #AP82156-1 50ug #AP82156-2 100ug #AP82156-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C chemokine receptor type 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CXC-R4,CD184
Accession No.	Uniprot ID:Q764M9
Uniprot	Q764M9
GeneID	396659;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CXCR4

Note: This product is for in vitro research use only