

C-X-C motif chemokine 15

Catalog No: #AP82157



Package Size: #AP82157-1 50ug #AP82157-2 100ug #AP82157-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 15
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lungkine, Small-inducible cytokine B15, Scyb15
Accession No.	Uniprot ID:Q9WVL7
Uniprot	Q9WVL7
GeneID	20309;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cxcl15

Note: This product is for in vitro research use only