

# Cytochrome P450 2E1

Catalog No: #AP82164



Package Size: #AP82164-1 50ug #AP82164-2 100ug #AP82164-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Cytochrome P450 2E1                                    |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Dog  |
| Other Names        | 4-nitrophenol 2-hydroxylase,CYP1IE1                    |
| Accession No.      | Uniprot ID:Q9MZY0                                      |
| Uniprot            | Q9MZY0   |
| GeneID             | 415128;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CYP2E1

Note: This product is for in vitro research use only