

N(G),N(G)-dimethylarginine dimethylaminohydrolase 2

Catalog No: #AP82167



Package Size: #AP82167-1 50ug #AP82167-2 100ug #AP82167-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N(G),N(G)-dimethylarginine dimethylaminohydrolase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DDAH-2,DDAHII,Dimethylargininase-2
Accession No.	Uniprot ID:Q99LD8
Uniprot	Q99LD8
GeneID	51793;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ddah2

Note: This product is for in vitro research use only