

DNA damage-inducible transcript 3 protein

Catalog No: #AP82168



Package Size: #AP82168-1 50ug #AP82168-2 100ug #AP82168-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA damage-inducible transcript 3 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DDIT-3,C/EBP zeta,C/EBP-homologous protein,CHOP,C/EBP-homologous protein 10,CHOP-10,CCAAT/enhancer-binding protein homologous protein,Growth arrest and DNA-damage-inducible protein GADD153,Chop,Chop10,Gadd153
Accession No.	Uniprot ID:Q62857
Uniprot	Q62857
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ddit3

Note: This product is for in vitro research use only