

25-hydroxycholesterol 7-alpha-hydroxylase

Catalog No: #AP82175



Package Size: #AP82175-1 50ug #AP82175-2 100ug #AP82175-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	25-hydroxycholesterol 7-alpha-hydroxylase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cytochrome P450 7B1, Hippocampal transcript 1 protein, HCT-1, Oxysterol 7-alpha-hydroxylase
Accession No.	Uniprot ID:Q60991
Uniprot	Q60991
GeneID	13123;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cyp7b1

Note: This product is for in vitro research use only