

cAMP-responsive element modulator

Catalog No: #AP82176



Package Size: #AP82176-1 50ug #AP82176-2 100ug #AP82176-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	cAMP-responsive element modulator
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q03061
Uniprot	Q03061
GeneID	25620;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Crem

Note: This product is for in vitro research use only