

CUB and zona pellucida-like domain-containing protein 1

Catalog No: #AP82180

Package Size: #AP82180-1 50ug #AP82180-2 100ug #AP82180-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CUB and zona pellucida-like domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CUB and ZP domain-containing protein 1,Estrogen-regulated protein 1,Uterus/ovary-specific protein 44,Erg1
Accession No.	Uniprot ID:Q9QZT0
Uniprot	Q9QZT0
GeneID	117179;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cuzd1

Note: This product is for in vitro research use only