

C-X-C chemokine receptor type 3

Catalog No: #AP82181



Package Size: #AP82181-1 50ug #AP82181-2 100ug #AP82181-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | C-X-C chemokine receptor type 3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | CXC-R3, Interferon-inducible protein 10 receptor, IP-10 receptor, CD183 |
| Accession No. | Uniprot ID:Q5MD61 |
| Uniprot | Q5MD61 |
| GeneID | 497018; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CXCR3

Note: This product is for in vitro research use only