

Cardiolipin synthase (CMP-forming)

Catalog No: #AP82187



Package Size: #AP82187-1 50ug #AP82187-2 100ug #AP82187-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cardiolipin synthase (CMP-forming)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CLS,Protein GCD10 homolog,C20orf155,CLS1
Accession No.	Uniprot ID:Q9UJA2
Uniprot	Q9UJA2
GeneID	54675;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CRLS1

Note: This product is for in vitro research use only