

Complement regulatory protein Crry

Catalog No: #AP82188



Package Size: #AP82188-1 50ug #AP82188-2 100ug #AP82188-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement regulatory protein Crry
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Antigen 5I2
Accession No.	Uniprot ID:Q63135
Uniprot	Q63135
GeneID	54243;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Crry

Note: This product is for in vitro research use only