

Cutaneous T-cell lymphoma-associated antigen 1

Catalog No: #AP82189



Package Size: #AP82189-1 50ug #AP82189-2 100ug #AP82189-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cutaneous T-cell lymphoma-associated antigen 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein cTAGE-1,Cancer/testis antigen 21.1,CT21.1,CTAGE2
Accession No.	Uniprot ID:Q9HC47
Uniprot	Q9HC47
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CTAGE1

Note: This product is for in vitro research use only