

Coenzyme Q-binding protein COQ10 homolog B, mitochondrial



Catalog No: #AP82192

Orders: order@signalwayantibody.com

Package Size: #AP82192-1 50ug #AP82192-2 100ug #AP82192-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Coenzyme Q-binding protein COQ10 homolog B, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q5I0I9
Uniprot	Q5I0I9
GeneID	301416;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Coq10b

Note: This product is for in vitro research use only