

Corticoliberin

Catalog No: #AP82195



Package Size: #AP82195-1 50ug #AP82195-2 100ug #AP82195-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Corticoliberin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Corticotropin-releasing factor,CRF,Corticotropin-releasing hormone
Accession No.	Uniprot ID:Q95MI6
Uniprot	Q95MI6
GeneID	280755;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CRH

Note: This product is for in vitro research use only