

Cystatin-B

Catalog No: #AP82197



Package Size: #AP82197-1 50ug #AP82197-2 100ug #AP82197-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Cystatin-B |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | CPI-B,Liver thiol proteinase inhibitor,Stefin-B,CST6,STFB |
| Accession No. | Uniprot ID:P04080 |
| Uniprot | P04080 |
| GeneID | 1476; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CSTB

Note: This product is for in vitro research use only