

Beta-defensin 103

Catalog No: #AP82198



Package Size: #AP82198-1 50ug #AP82198-2 100ug #AP82198-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-defensin 103
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Beta-defensin 3,BD-3,Defensin, beta 103,Defensin-like protein,BD3,DEFB103,DEFB3
Accession No.	Uniprot ID:P81534
Uniprot	P81534
GeneID	414325;55894;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DEFB103A

Note: This product is for in vitro research use only