

Coenzyme Q-binding protein COQ10 homolog B, mitochondrial



Catalog No: #AP82200

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AP82200-1 50ug #AP82200-2 100ug #AP82200-3 1mg

Description

Product Name	Coenzyme Q-binding protein COQ10 homolog B, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q9H8M1
Uniprot	Q9H8M1
GeneID	80219;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: COQ10B

Note: This product is for in vitro research use only