

Cullin-5

Catalog No: #AP82204



Package Size: #AP82204-1 50ug #AP82204-2 100ug #AP82204-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cullin-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CUL-5
Accession No.	Uniprot ID:Q9D5V5
Uniprot	Q9D5V5
GeneID	75717;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cul5

Note: This product is for in vitro research use only