

# Gap junction delta-2 protein

Catalog No: #AP82205

Package Size: #AP82205-1 50ug #AP82205-2 100ug #AP82205-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gap junction delta-2 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Connexin-36,Cx36,Gap junction alpha-9 protein,GJA9
Accession No.	Uniprot ID:Q866T7
Uniprot	Q866T7
GeneID	282577;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:GJD2

Note: This product is for in vitro research use only